

# Motor X

## Smart Implant Torque Wrench

One-handed operation

Rapid treatment

Stable and accurate



# Rapid treatment

More than 2 times faster than using manual torque wrench

Motor X





VS

Traditional torque wrench



## Adjustable direction of rotation



Motor X   
 40 rpm  
40 N·cm

Forward to implant

Motor X   
 20 rpm  
20 N·cm

Reverse to remove


# Security intelligence

Automatically stop when the set torque is encountered under the forward rotation mode



# High accuracy

The error of torque and speed is controlled within 10%



## LCD digital display

Show every function including torque, speed and direction of rotation



# Long battery life

120 minutes of continuous operation time  
when fully charged



1900mAh

# Packing list:

## 16 types of screwdrivers





# One-handed operation

Provides a wider view of operation site



# Ergonomic design

Make clinicians to operate easily

